# INHEALTH

TECHNOLOGIES®

## Blom-Singer® Adjustable Tracheostoma Valve II

### INSTRUCTIONS FOR USE

**Blom-Singer®**

voice restoration systems

37433-01H Effective January 2020 / 2020 年 1 月生效 / 2020年1月発効 / Në fuqi nga janari 2020 / 2020 اعتبارًا من شهر يناير / В сила от януари 2020 г. / Stupa na snagu u siječnju 2020. / V platnosti od ledna 2020 / Gældende fra januar 2020 / Met ingang van januari 2020 / Kehtib alates jaanuarist 2020 / Voimassa tammikuusta 2020 alkaen / Date d'entrée en vigueur : janvier 2020 / Gültig ab Januar 2020 / Σε ισχύ από τον Ιανουάριο του 2020 / 2020 בתוקף החל מינואר / Hatályos 2020. januártól / Gildir frá janúar 2020 / Valido da gennaio 2020 / 2020년 1월 발효 / Spēkā no 2020. gada janvāra / Įsigalioja 2020 m. sausio mėn. / Gjelder fra januar 2020 / Wchodzi w życie w styczniu 2020 r. / Vigente em janeiro de 2020 / Cu intrare în vigoare din ianuarie 2020 / Версия от января 2020 г. / Važi od januara 2020. / Účinné od januára 2020 / Efectivo a partir de enero de 2020 / Gäller från januari 2020 / Geçerlilik Tarihi: Ocak 2020

Blom-Singer® and InHealth Technologies® are registered trademarks of Freudenberg Medical, LLC.

Patent(s): [www.inhealthpatents.com](http://www.inhealthpatents.com)




Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.




## INSTRUCTIONS FOR USE




使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

## DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagrame / Схеми / Crteži / Grafické znázornenia / Diagramas / Diagram / Diyagramlar

## TABLE OF CONTENTS

目录 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /  
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières /  
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /  
 Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins /  
 Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Adjustable Tracheostoma Valve II	7	HU	Állítható tracheostoma szelep II	64
ZH	可调节式气管造口阀 II	10	IS	Stillanlegur ventill fyrir barkarauf II	67
JA	調節可能気管孔バルブ II	13	IT	Valvola regolabile per tracheostomia II	70
SQ	Valvula e rregullueshme e trakeostomës II	17	KO	조절 가능 기관절개창 밸브 II	74
20	صمام الفقرة الرغامية II القابل للتعديل	AR	LV	Pielāgojamais traheostomas vārsts II	78
BG	Регулируем трахеостомен клапан II	24	LT	Reguliuojamas traheostomos vožtuvas II	81
HR	Prilagodljivi traheostomni ventil II	28	NO	Justerbar tracheostoma ventil II	85
CS	Nastavitelný tracheostomický ventil II	31	PL	Regulowana zastawka tracheostomijna II	88
DA	Justerbar tracheostomiventil II	35	PT	Válvula ajustável para traqueostoma II	92
NL	Aanpasbare tracheostomaklep II	38	RO	Valvă reglabilă pentru traheostomă II	96
ET	Reguleeritav trahheestoomiklapp II	42	RU	Регулируемый клапан трахеостомы II	99
FI	Säädettävä tracheostoomaläppä II	45	SR	Podesiv traheostomni ventil II	103
FR	Clapet pour trachéostomie réglable II	48	SK	Nastavitelný tracheostomický ventil II	107
DE	Anpassbares Tracheostomaventil II	52	ES	Válvula para traqueostoma ajustable II	110
EL	Ρυθμιζόμενη βαλβίδα τραχειοστομίας II	56	SV	Justerbar tracheostomiventil II	114
60	II שתום פיום קנה מתכוון	HE	TR	Ayarlanabilir Trakeostoma Valfi II	118

## ENGLISH

**BLOM-SINGER® ADJUSTABLE TRACHEOSTOMA VALVE II****PRODUCT DESCRIPTION**

The Blom Singer Adjustable Tracheostoma Valve II (ATSV II) is a pre-assembled hands-free valve, reusable for single-patient only (diagram 1). The device consists of removable diaphragm/faceplate (1A) and a body (1B). The assembled two-piece valve inserts into a tracheostoma housing (diagrams 6 and 7) that acts as a holder for the assembled device. The ATSV II consists of the following components:

- **Diaphragm/Faceplate:** The faceplate has a groove that aligns with the peg on the body (diagram 3a) to hold these two components together. A replacement Diaphragm/Faceplate is available separately.
- **Diaphragm Flap Valve:** The silicone diaphragm flap valve is permanently attached to the bottom of the faceplate (diagram 3b). It is designed to facilitate occlusion of the stoma to enable hands-free speech with a voice prosthesis.
- **Crossbar:** The crossbar (diagram 1a) is designed to prevent clothing from blocking the faceplate of the device and for holding the optional HumidiFilter system (supplied separately) in place. It permits adjustment of the diaphragm/flap valve. The crossbar should not be used as a handle for inserting or removing the ATSV II from the housing. Do not cut or remove the crossbar, or the device will become unusable.
- **Body:** The body is designed to hold the diaphragm/faceplate. It has a peg (diagram 2a) that aligns with the groove in the diaphragm/faceplate.

A Blom-Singer Adhesive or Tracheostoma Housing (supplied separately) is necessary for attachment of the ATSV II over the tracheostoma. For information about how to use the housings, please refer to the instructions for use provided with the selected housing. The ATSV II is compatible with the Blom-Singer® HumidiFilter® System. The ATSV II HumidiFilter Cap fits over the device.

**INDICATIONS** (Reasons to prescribe the device or procedure)

The Blom-Singer Adjustable Tracheostoma Valve II (ATSV II) is a medical device intended for use only by laryngectomees who have undergone surgical voice restoration and use a voice prosthesis. When used in conjunction with a voice prosthesis, the ATSV II eliminates the need for digital stoma occlusion to produce voice.

**CONTRAINDICATIONS** (Reasons that make it inadvisable to prescribe the particular device or procedure)

The ATSV II should not be placed over a tracheostomal airway if its presence obstructs respiration.

**WARNINGS AND PRECAUTIONS**

Freudenberg Medical relies on the qualified, trained medical professional to instruct patients on the use and care of this device and to provide patients with the instructions for use supplied with the product. Not following these instructions could lead to skin or lung irritation. If there are tears, cracks, or structural damage to the device, discontinue use. Extreme care should be used to prevent the device from entering the tracheostoma.

Breathing

Breathing through the ATSV II device should be comfortable and easy. If breathing is difficult, remove the ATSV II device and consult your medical professional for evaluation. When you begin using the ATSV II device, a slight increase in airflow resistance may be felt. This resistance is a normal

function of the ATSV II device that may require a period of time to adjust. The ATSV II may require temporary removal during instances of extreme physical activity or heavy breathing.

### Coughing

The ATSV II and HumidiFilter Cap (if applicable) should be removed from the tracheostoma housing prior to coughing.

### Cleaning

Do not use solvents or petroleum-based products for cleaning or lubricating the device. These materials may damage the device or cause it to not work properly. Never use facial or toilet tissue to dry the device as particles of lint or fabric may be deposited onto the device that can be inhaled into the airway upon reinsertion. Either air-dry or use a lint-free material to dry the device. Never clean the ATSV II, or insert any cleaning device, while the ATSV II is in position in the neck. Only clean the device after it has been removed.

### Sleeping

The ATSV II must be removed before sleeping to avoid the potential for airway obstruction.

## **INSTRUCTIONS FOR USE**

*The following procedural instructions are provided by Eric D. Blom, Ph.D. Please refer to the diagrams located at the front of this instruction manual.*

Your clinician will instruct you on the routine use and care of the device. The ATSV II, device components, tracheostoma housing, and user's hands should be clean to prevent contaminants from entering into the body. When clinicians are handling, gloved hands are recommended. Removal and attachment of any devices worn over the tracheostoma should only be done with the user positioned in front of a mirror with a bright light focused directly on the stoma. Before inserting an ATSV II device, the tracheostoma housing must first be attached to the skin surrounding the tracheostoma. Follow the attachment steps described in the tracheostoma housing Instructions For Use provided with the selected product.

### Insertion

To insert the ATSV II, grasp its outer rim and insert it partially into the lower rim of the previously applied tracheostoma housing (diagram 8). Gently push the ATSV II until it snaps into place. A properly inserted ATSV II will fit securely inside the housing (diagram 9). Do not use the crossbar as a handle to insert the ATSV II.

### Removal

The ATSV II is designed for removal without requiring removal of the tracheostoma housing from the skin. To remove the device, hold down the previously applied tracheostoma housing with one hand while grasping the outer rim of the device with your other hand. Pull the device carefully until it releases from the housing (diagram 10). Do not use the crossbar as a handle to remove the ATSV II.

### Cleaning

The ATSV II should be removed and soaked in warm water for 2-3 minutes to help break down any mucus that may have accumulated in the device. Carefully clean the surfaces of the diaphragm/faceplate and body with a small brush or cotton tip applicator and mild soap. Do not use strong detergents or any other chemical that may be damaging to the device. Failure to clean the device thoroughly may interfere with the proper function of the diaphragm flap valve. Rinse with water and either air-dry or use a lint-free material to dry the device. Never use facial or toilet tissue to dry the ATSV II as particles of lint or fabric may be deposited onto the device that can be inhaled into the airway upon reinsertion. Do not clean or place the ATSV II in water while the HumidiFilter Cap is attached.

### ATSV II Function

When properly adjusted, the diaphragm flap valve remains in an open position during normal, quiet breathing and routine physical activity. As the user exhales to speak, the airflow pushes the valve to the closed position to provide a seal during speech. When the user stops talking or inhales, the diaphragm flap valve reopens and returns to its resting position. With the extremely high pressure of a cough, the diaphragm flap valve may be forced partially through the opening in the faceplate. This may be accompanied a “snapping” sound. If this happens, remove the ATSV II from the tracheostoma housing and reposition the diaphragm flap valve with a clean fingertip.

### ATSV II Sensitivity Adjustments

The ATSV II is designed to accommodate the speech and breathing requirements of each individual user. Valve sensitivity adjustment is most easily accomplished by rotating the faceplate when the ATSV II is removed from the tracheostoma housing. Clockwise rotation adjusts the diaphragm flap valve to a more closed position to make it close easily during speech (diagram 11). Counterclockwise rotation adjusts the diaphragm flap valve to a more open position thus requiring greater effort (respiratory airflow) to close it to speak (diagram 12). Most users set the ATSV II to a near-closed position during speech and adjust it to a wide-open position during increased respiratory exertion (e.g., during physical exercise). Intermediate positions can be selected in accordance with varying speech and breathing requirements of the user. To adjust the ATSV II while it is in position in the adhesive housing, the user should stabilize the base of the ATSV II between the index finger and thumb of one hand while simultaneously rotating the faceplate with the other hand. Your clinician/physician can assist you in determining the appropriate valve sensitivity settings for you.

### **Diaphragm/Faceplate Replacement**

If the diaphragm flap valve becomes worn or damaged, it should be replaced. A replacement diaphragm/faceplate is available separately from InHealth Technologies. Always remove the ATSV II from the tracheostoma housing before disassembling it or replacing the diaphragm/faceplate. To replace the diaphragm/faceplate from the ATSV II body, insert a fingernail or other small, flat object into the slot located on the side of the ATSV II. Gently twist to detach the diaphragm/faceplate from the body. Discard the used diaphragm/faceplate. Align the peg on the body (diagram 13) with the groove in the rim of the new diaphragm/faceplate replacement and “snap” these two components together (diagram 14). The diaphragm flap valve should curl over the small “plastic ledge” (diagram 13a) inside the ATSV II. Confirm proper assembly of the ATSV II and the diaphragm flap valve prior to use.

### **TROUBLESHOOTING**

**Problem:** Air leak through or around diaphragm flap valve during speech.

**Solution:** Disassemble ATSV II and thoroughly clean any mucus that may prevent the diaphragm flap valve from closing during speech. Replace diaphragm/faceplate if it has become worn or does not work properly.

**Problem:** Diaphragm flap valve is not easily closed during speech.

**Solution:** Rotate faceplate clockwise to adjust diaphragm flap valve to a more closed position.

**Problem:** Diaphragm/faceplate closes inadvertently.

**Solution:** Rotate faceplate counterclockwise to adjust the diaphragm flap valve setting to a more open position.

If you experience difficulty using the ATSV II device, discontinue use and consult your clinician for evaluation.

## SPECIAL STORAGE AND/OR HANDLING CONDITIONS

There are no special storage and/or handling conditions for this device.

## DISPOSAL INSTRUCTIONS

This device is not biodegradable and can be contaminated when used. Carefully dispose of the device per local guidelines.

## ORDERING INFORMATION

### USA

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: [order@inhealth.com](mailto:order@inhealth.com) ORDER ONLINE: [www.inhealth.com](http://www.inhealth.com) POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

### INTERNATIONAL

Please contact our customer service department for a distributor referral.

## PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact [productcomplaints@inhealth.com](mailto:productcomplaints@inhealth.com).

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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**SYMBOLS GLOSSARY (USA only)**

<b>SYMBOL</b>	<b>STANDARD No.</b>	<b>SYMBOL TITLE &amp; REF No.</b>	<b>MEANING OF SYMBOL</b>
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Batch code; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the packaged has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	ISO 7000:2019	Single patient multiple use; 3706	To indicate that the medical device may be used multiple times (multiple procedures) on a single patient

# INHEALTH

## TECHNOLOGIES®

---

### UNITED STATES

1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
Toll-Free: 800.477.5969  
Toll-Free Fax: 888.371.1530

### INTERNATIONAL

Tel: +1.805.684.9337  
Fax: +1.805.684.8594  
E-mail: [order@inhealth.com](mailto:order@inhealth.com)

### ORDER ONLINE

[www.inhealth.com](http://www.inhealth.com)

**Freudenberg Medical, LLC**  
1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
[www.inhealth.com](http://www.inhealth.com)

**EMERGO EUROPE**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands